



Reed
Biotech



FAST STEP ELISA KIT

FastStep™
One-Step Method

Sample Volume

Assay Time

Washing Times

Incubation Times

Operation Difficulty

50μL

1.5h

5 times

1 time

Easy

VS

Sample Volume

Assay Time

Washing Times

Incubation Times

Operation Difficulty

100μL

3.5h

9 times

3 times

Relatively difficult

Traditional
Three-Step
Method


- Sample Saving

Requires only 50uL sample–saves 50% of the sample volume, compared with the traditional three-step method.
- Time Saving

Reaction time is only 1.5h, 2hs shorter than traditional three-step method.
- Simpler Operation

Only a single incubation is needed, significantly improving experimental reproducibility.
- Higher Efficiency

Greatly simplifies experimental procedures to enhance efficiency.



FastStep™ HOT TARGETS				
Cat. NO	Product Name	Reaction time	Sensitivity	Detection range
FS10109	FS-Cor(Cortisol) ELISA Kit	0.75H	0.17 ng/mL	0.31-20ng/mL
FS10118	FS-Pg(Progesterone) ELISA Kit	0.75H	42.97 pg/mL	78.13-5000 pg/mL
FS10133	FS-ALD(Aldosterone) ELISA Kit	0.75H	14.61 pg/mL	31.25-2000pg/mL
FS10177	FS-E2(Estradiol) ELISA Kit	0.75H	11.25 pg/mL	12.5-800pg/mL
FS1033H	FS-Human EGF(Epidermal Growth Factor) ELISA Kit	1.5H	0.09 pg/mL	1.56-100pg/mL
FS1049R	FS-Rat IL-1R2(Interleukin 1 Receptor Type II) ELISA Kit	1.5H	5.57 pg/mL	15.63-1000pg/mL
FS1094H	FS-Human TIMP-1(Tissue Inhibitors of Metalloproteinase 1) ELISA Kit	1.5H	31.25 pg/mL	78.13-5000pg/mL
FS1094M	FS-Mouse TIMP-1(Tissue Inhibitors of Metalloproteinase 1) ELISA Kit	1.5H	25.41 pg/mL	78.13-5000pg/mL
FS1100M	FS-Mouse IL-6R(Interleukin 6 Receptor) ELISA Kit	1.5H	20.61 pg/mL	39.06-2500pg/mL
FS1305H	FS-Human HBP/AZU1(Azurocidin 1) ELISA Kit	1.5H	22.30 pg/mL	62.5-4000pg/mL
FS1577H	FS-Human IL-2Rα(Interleukin 2 Receptor Alpha) ELISA Kit	1.5H	1.26 pg/mL	15.63-1000pg/mL
FS1739H	FS-Human SCF(Stem Cell Factor) ELISA Kit	1.5H	3.76 pg/mL	31.25-2000pg/mL
FS1828H	FS-Human FABP3(Fatty Acid Binding Protein 3, Muscle and Heart) ELISA Kit	1.5H	0.10 ng/mL	0.16-10ng/mL
FS1917M	FS-Mouse PGRN(progranulin) ELISA Kit	1.5H	0.04 ng/mL	0.16-10ng/mL
FS2005H	FS-Human SCCA1(Squamous Cell Carcinoma Antigen 1) ELISA Kit	1.5H	10.01 pg/mL	78.13-5000pg/mL
FS2103H	FS-Human NCAM/CD56(Neural Cell Adhesion Molecule) ELISA Kit	1.5H	0.07 ng/mL	0.16-10ng/mL
FS2202H	FS-Human PGA/Pg I (Pepsinogen A) ELISA Kit	1.5H	0.23 ng/mL	0.39-25ng/mL
FS2203H	FS-Human PGC/Pg II (Pepsinogen C) ELISA Kit	1.5H	24.79 pg/mL	31.25-2000pg/mL
FS2216H	FS-Human CYFRA21-1(Cytokeratin Fragment Antigen 21-1) ELISA Kit	1.5H	0.15 ng/mL	0.625-40ng/mL
FS2731H	FS-Human MCP-1(Monocyte Chemotactic Protein 1) ELISA Kit	1.5H	3.48 pg/mL	15.63-1000 pg/mL
FS2739M	FS-Mouse LEP(Leptin) ELISA Kit	1.5H	29.92 pg/mL	78.13-5000pg/mL
FS2742H	FS-Human C-P(C-Peptide) ELISA Kit	1.5H	18.24 pg/mL	31.25-2000pg/mL
FS2746H	FS-Human MYO(Myoglobin) ELISA Kit	1.5H	0.94 ng/mL	1.56-100ng/mL
FS2796M	FS-Mouse MMP-9(Matrix Metalloproteinase 9) ELISA Kit	1.5H	0.03 ng/mL	0.31-20ng/mL
FS2796R	FS-Rat MMP-9(Matrix Metalloproteinase 9) ELISA Kit	1.5H	0.05 ng/mL	0.31-20ng/mL
FS2831H	FS-Human ICAM-1/CD54(intercellular adhesion molecule 1) ELISA Kit	1.5H	18.75 pg/mL	31.25-2000pg/mL
FS2831R	FS-Rat ICAM-1/CD54(Inter cellular Adhesion Molecule 1) ELISA Kit	1.5H	0.01 ng/mL	0.31-20ng/mL
FS2837H	FS-Human ADP/Acrp30(Adiponectin) ELISA Kit	1.5H	6.94 pg/mL	62.5-4000pg/mL
FS2840H	FS-Human NGAL(Neutrophil Gelatinase Associated Lipocalin) ELISA Kit	1.5H	0.16 ng/mL	0.31-20ng/mL
FS2854M	FS-Mouse CRP(C-ReactiveProtein) ELISA Kit	1.5H	0.38 ng/mL	0.63-40ng/mL
FS2872H	FS-Human FE(Ferritin) ELISA Kit	1.5H	0.04 ng/mL	0.31-20ng/mL
FS2873H	FS-Human IgG(Immunoglobulin G) ELISA Kit	1.5H	1.45 ng/mL	3.13-200ng/mL
FS2873M	FS-Mouse IgG(Immunoglobulin G) ELISA Kit	1.5H	0.09 ng/mL	1.56-100ng/mL
FS2892H	FS-Human Ang- I (Angiotensin I) ELISA Kit	0.75H	0.12 ng/mL	0.63-40ng/mL
FS2908H	FS-Human IGFBP-1(insulin-like growth factors binding protein 1) ELISA Kit	1.5H	0.09 ng/mL	0.16-10ng/mL
FS2970R	FS-Rat C1qR1(Complement Component 1q Receptor 1) ELISA Kit	1.5H	0.09 ng/mL	0.16-10ng/mL
FS3008H	FS-Human NSE ELISA Kit	1.5H	15.63 pg/mL	78.13-5000 pg/mL
FS3064H	FS-Human PDGF-BB(Platelet Derived Growth Factor BB) ELISA Kit	1.5H	25.92 pg/mL	31.25-2000pg/mL
FS3122R	FS-Rat BMG/β2-MG(Beta-2-Microglobulin) ELISA Kit	1.5H	3.14 pg/mL	31.25-2000pg/mL
FS3130H	FS-Human CA242(Pancreatic Carcinoma Markers-CA242) ELISA Kit	1.5H	0.05 IU/mL	0.16-10IU/mL
FS3153H	FS-Human INS(Insulin) ELISA Kit	1.5H	0.16 μIU/mL	0.31-20μIU/mL
FS3186H	FS-Human IL-6(Interleukin 6) ELISA Kit	1.5H	0.81 pg/mL	1.56-100pg/mL
FS3197H	FS-Human ENG(Endoglin) ELISA Kit	1.5H	0.03 ng/mL	0.16-10ng/mL
FS3393R	FS-Rat sCD86(Soluble Cluster of Differentiation 86) ELISA Kit	1.5H	8.93 pg/mL	15.63-1000pg/mL
FS10113	FS-Cort(Corticosterone) ELISA Kit	0.75H	2.76 ng/mL	3.13-200ng/mL