

Cystatin C Mouse Monoclonal Antibody(3B12)

Catalog No: RA10118

Basic Information

Host species	Mouse
Applications	WB,ELISA
Species Cross-Reactivity	H
Specificity	Antibody can detects endogenous Cystatin C protein.
Recommended dilutions	WB: 1:1,000-2,000 ELISA: 1:10,000-20,000 Optimal dilutions should be determined by the end user.

Applications

Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1 mg/ml
Clonality	Monoclonal

Background

AD8, Cystatin3, Gamma trace, HCCAA

Alternative Names	_x0007_
	x0007
	x0007
	x0007
	x0007
	x0007
	x0007
	x0007
	x0007
	x0007
Observed band	14

Human Gene ID	1471
Human Swiss-Prot Number	P01034
Background	