

alpha-Syneuclein Mouse Monoclonal Antibody(2H5)

Catalog No: RA10278

Basic Information

Host species	Mouse
Applications	ELISA
Species Cross-Reactivity	H
Specificity	The antibody detects endogenous alpha-Syneuclein proteins.
Recommended dilutions	ELISA 1:2,000-10,000 Optimal dilutions should be determined by the end user.

Applications

Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1 mg/ml
Clonality	Monoclonal

Background

Alternative Names	PARK1, SYN, SNCA, MGC105443,NACP
Observed band	14
Human Gene ID	6622
Human Swiss-Prot Number	P37840
Background	Genetic alterations of SNCA resulting in aberrant polymerization into fibrils, are associated with several neurodegenerative diseases (synucleinopathies). SNCA fibrillar aggregates represent the major non A-beta component of Alzheimer disease amyloid plaque, and a major component of Lewy body inclusions.