

# 14-3-3 eta Mouse Monoclonal Antibody(8H2)

## Catalog No: RA10283

### Basic Information

Host species	Mouse
Applications	ELISA
Species Cross-Reactivity	H, M, R
Specificity	The antibody detects endogenous 14-3-3 eta proteins.
Recommended dilutions	ELISA 1:5,000-10,000 Optimal dilutions should be determined by the end user.

### Applications

Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1 mg/ml
Clonality	Monoclonal

### Background

Alternative Names	YWHAH; YWHA1; 14-3-3 protein eta; Protein AS1
Observed band	28
Human Gene ID	7533
Human Swiss-Prot Number	Q04917
Background	tyrosine 3-monooxygenase/tryptophan 5-monooxygenase activation protein eta(YWHAH) Homo sapiens This gene product belongs to the 14-3-3 family of proteins which mediate signal transduction by binding to phosphoserine-containing proteins. This highly conserved protein family is found in both plants and mammals, and this protein is 99% identical to the mouse, rat and bovine orthologs. This gene contains a 7 bp repeat sequence in its 5' UTR, and changes in the number of this repeat have been associated with early-onset schizophrenia and psychotic bipolar disorder.