

# MFAP4 Mouse Monoclonal Antibody(2C12)

## Catalog No: RA10286

### Basic Information

Host species	Mouse
Applications	ELISA, IHC
Species Cross-Reactivity	H, M, R
Specificity	The antibody detects endogenous MFAP4 proteins.
Recommended dilutions	ELISA 1:5,000-10,000 IHC 1:100-200 Optimal dilutions should be determined by the end user.

### Applications

Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1 mg/ml
Clonality	Monoclonal

### Background

Alternative Names	Microfibril associated glycoprotein 4
Observed band	
Human Gene ID	4239
Human Swiss-Prot Number	P55083
Background	MFAP4 is deleted in the Smith-Magenis syndrome.