

Ubiquitin Mouse Monoclonal Antibody(2G3)

Catalog No: RA10290

Basic Information

Host species	Mouse
Applications	WB,IHC
Species Cross-Reactivity	ALL
Specificity	Antibody can detects endogenous Ubiquitin protein.
Recommended dilutions	WB: 1:1,000-2,000 IHC: 1:200-500 Optimal dilutions should be determined by the end user.


Applications

Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1 mg/ml
Clonality	Monoclonal

Background

Alternative Names	UBB, UBC,UBA80,UBCEP1,CEP52
Observed band	ALL
Human Gene ID	7314
Human Swiss-Prot Number	PAN
Background	Ubiquitin is a conserved polypeptide unit that plays an important role in the ubiquitin-proteasome pathway. Ubiquitin can be covalently linked to many cellular proteins by the ubiquitination process, which targets proteins for degradation by the 26S proteasome.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma Tissue using Ubiquitin Mouse Monoclonal antibody diluted at 1:500.