

# Histone H3(tri methyl K79) Mouse Monoclonal Antibody(3G3)

## Catalog No: RA10369

### Basic Information

Host species	Mouse
Applications	WB, IP, IHC
Species Cross-Reactivity	H, R, M
Specificity	Antibody can detects endogenous Histone H3 (tri methyl K79) protein.
Recommended dilutions	WB: 1:500-2,000      IP:1:200      IHC: 1:200-500 Optimal dilutions should be determined by the end user.


### Applications

Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1 mg/ml
Clonality	Monoclonal


### Background

Alternative Names	H3 histone antibody, HIST1H3A antibody, Histone cluster 1, H3a antibody
Observed band	15
Human Gene ID	8290
Human Swiss-Prot Number	P68431
Background	Histone H3 is one of the five main histone proteins involved in the structure of chromatin in eukaryotic cells. Core component of nucleosome. Nucleosomes wrap and compact DNA into chromatin, limiting DNA accessibility to the cellular machineries which require DNA as a template. Histones thereby play a central role in transcription regulation, DNA repair, DNA replication and chromosomal stability.

## Selected Validation Data



Western blot analysis of 1) HeLa, 2) Rat Testis tissue, 3) Raw264.7 with Histone H3 (tri methyl K79) (3G3) Mouse MAb diluted at 1:2000.



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma using Histone H3(tri methyl K79) Mouse mAb diluted at 1:500.