

Zip9 Mouse mAb(Mix-mA)

Catalog No: RA10411

Basic Information

Host species	Mouse
Applications	WB,IHC
Species Cross-Reactivity	H,M,R
Specificity	Antibody can detects recombinant Zip9 protein.
Recommended dilutions	WB:1:1,000-2,000 IHC:1:200-500 Optimal dilutions should be determined by the end user.

Applications

Formulation	Antigen Affinity Purified IgG in PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1 mg/ml
Clonality	Monoclonal

Background

Alternative Names	SLC39A9, Zinc transporter zip9, solute carrier family 39 member 9
Observed band	33
Human Gene ID	55334
Human Swiss-Prot Number	Q9NUM3
Background	function:May act as a zinc-influx transporter.,similarity:Belongs to the ZIP transporter (TC 2.A.5) family.