

# EBFP Mouse Monoclonal Antibody(Mix-mA™)

## Catalog No: RA10434

### Basic Information

Host species	Mouse
Applications	WB, IHC, IP
Species Cross-Reactivity	N/A
Specificity	The EBFP Mouse Polyclonal Antibody can detects EBFP & EBFP Tag protein.
Recommended dilutions	WB: 1:5,000    IHC: 1:200    IP: 1:200 Optimal dilutions should be determined by the end user.

### Applications

Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1 mg/ml
Clonality	Monoclonal

### Background

Alternative Names	BFP, EBFP, Blue Fluorescent Protein, Enhanced Blue Fluorescent Protein
Observed band	N/A
Human Gene ID	
Human Swiss-Prot Number	N/A
Background	Enhanced Blue Fluorescent Protein (EBFP) emits a strong blue fluorescence. EBFP can be used as protein "tags" to study the sub-cellular localization of proteins and/or their translocation upon stimulation or as markers for transfection in transient and stable expression systems.

Selected Validation Data



Western blot analysis of recombinant EBFP protein with EBFP  
Mouse Monoclonal Antibody(Mix-mA™) diluted at 1:5,000.