

VEGF Rabbit Polyclonal Antibody(VEGF)

Catalog No: RA20266

Basic Information

Host species	Rabbit
Applications	IHC
Species Cross-Reactivity	H,R,M
Specificity	Antibody can detects endogenous VEGF protein.
Recommended dilutions	IHC: 1:200-500 Optimal dilutions should be determined by the end user.


Applications

Formulation	Antigen Affinity Purified IgG in PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1 mg/ml
Clonality	Polyclonal


Background

	VEGFA antibody, VPF antibody, VEGF120 antibody, Vascular endothelial growth factor antibody
Alternative Names	/ _x0007_ _x0007_ _x0007_Immunohistochemical analysis of paraffin-embedded Human Liver Carcinoma Tissue using VEGF (TDY649) Rabbit pAb diluted at 1:500. _x0007_Immunohistochemical analys
Observed band	23,45
Human Gene ID	7422
Human Swiss-Prot Number	P15692
Background	

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded Human Liver Carcinoma Tissue using VEGF Rabbit pAb diluted at 1:500.



Immunohistochemical analysis of paraffin-embedded Human Stomach Carcinoma Tissue using VEGF Rabbit pAb diluted at 1:500.