

# Luciferase Rabbit Polyclonal Antibody

## Catalog No: RA20321

### Basic Information

Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Firefly
Specificity	The Luciferase antibody can detect endogenous and recombinant Luciferase protein.
Recommended dilutions	WB: 1:1,000-2,000 Optimal dilutions should be determined by the end user.

### Applications

Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1 mg/ml
Clonality	Polyclonal

### Background

Alternative Names	
Observed band	
Human Gene ID	6737
Human Swiss-Prot Number	P19474
Background	

Selected Validation Data



Western blot analysis of Recombinant FireFly protein using  
Luciferase Rabbit pAb diluted at 1:2,000.