

mStrawberry Mouse Polyclonal Antibody

Catalog No: RA20360

Basic Information

Host species	Mouse
Applications	WB, IHC, IF, IP
Species Cross-Reactivity	N/A
Specificity	The mStrawberry Mouse Polyclonal Antibody can detect mStrawberry & mStrawberry Tag protein.
Recommended dilutions	WB: 1:5,000 IHC: 1:200 IF: 1:200 IP: 1:200 Optimal dilutions should be determined by the end user.

Applications

Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1 mg/ml
Clonality	Polyclonal

Background

Alternative Names	mStrawberry
Observed band	N/A
Human Gene ID	
Human Swiss-Prot Number	N/A
Background	mStrawberry is a bright, monomeric red fluorescent protein which was developed in Dr. Roger Tsien's lab by directed mutagenesis of mRFP.

Selected Validation Data



Western blot analysis of mStawberry recombinant protein with
mStrawberry Mouse Polyclonal Antibody diluted at 1:5,000.