

SV40 T-antigen Rabbit Polyclonal Antibody

Catalog No: RA20403

Basic Information

Host species	Rabbit
Applications	WB
Species Cross-Reactivity	
Specificity	Antibody can detects Recombinant SV40 T-antigen proteins.
Recommended dilutions	WB: 1:10,000-20,000 Optimal dilutions should be determined by the end user.

Applications

Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1 mg/ml
Clonality	Polyclonal

Background

Alternative Names	
Observed band	KD
Human Gene ID	0
Human Swiss-Prot Number	0
Background	.

Selected Validation Data



Western blot analysis of Recombinant protein using SV40
T-antigen Rabbit pAb diluted at 1:10,000.