

# E2-Tag Mouse Monoclonal Antibody(4F8)

## Catalog No: RTA31

### Basic Information

Host species	Mouse
Applications	WB
Species Cross-Reactivity	N/A
Specificity	The E2 tag antibody can recognize E2 tag fusion proteins.
Recommended dilutions	WB: 1:5,000 Optimal dilutions should be determined by the end user.

### Applications

Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1 mg/ml
Clonality	Monoclonal

### Background

Alternative Names	
Observed band	N/A
Human Gene ID	N/A
Human Swiss-Prot Number	N/A
Background	E2-tag peptide GVSSTSSDFRDR conjugated to KLH. E2- Tag antibody can recognize C-terminal, internal, and N-terminal E2-tagged proteins.

Selected Validation Data

